


Polyester Braided Tapes



Offered By TDA For

- **Umibilical tie off**
- **Retraction**
- **Ligation**

The braided polyester tapes are made from continuous filament high strength polyester yarns. The tapes have very good handling characteristics and smooth surfaces to facilitate their ease of use.

Properties	HTPETBR03S	HTPETBR05S
Width, mm	3	5
Thickness, mm	0.49	0.49
Strength, kg	29	55
Knot Strength, kg	16	28
Elongation, %	32	32
Yield, m/g	1.25	0.65

DISCLAIMER: The information supplied in this document is for guidance only and should not be construed as a warranty. All implied warranties are expressly disclaimed, including without limitation any warranty of merchantability and fitness for use. All users of the material are responsible for assuring that it is suitable for their needs, environment and end use. All data is subject to change as TDA deems appropriate.